

Blood Tubular Separator Series



CHINA

No.13093-06-02 Victory Community Baita District Liaoyang China

PHILIPPINES

Suite 310 Intramuros Corporate Plaza Recoletos St., Intramuros Manila Philippines

USA

6900 46th Street Kenosha WI 53144

KOREA

5-4, Songnim-Ro 48 beon-Gil yuseong-Gu, Daejon Korea



sales@crown-machinery.com



+86 186 4192 8887



China +86 419 2883215 Philippines +63 2 2549729 USA +1 262 656 7680 Korea +82 42 825 4260



+86 419 2122006

Base on 30 years for tubular separator design and market experience, Crown Machinery designed and manufacture the China made tubular separator, earn a reputation and good feedback from global customers.

Blood powder (animal) is the collection of quarantine health of fresh blood, through high temperature and pressure dispose, bleaching. low temperature drying and other processes made of high quality protein feed, with digestibility and safety, is high quality Animal protein source, suitable for livestock and aquatic animals, the daily food preparation. Centrifugal a membrane filter pre-treatment separator as equipment., which separate and extract the blood cell and plasma.

Crown-Machinery Tubular Separator



TSX SERIES SEPARATION TYPE

The Crown Machinery TSX series of tubular separator can be configured as a three phase solutions (Plasma and Blood cell). with the advantage of widely range in animal or human blood, accurate in plasma and blood cell separation, high yield of plasma, special design ensure the blood cell will not be destroied.

SWING TYPE

Considering of female workers , the swing type can easily pull out the bowl and fit it on the bowl wheel holder. Contrasting with the Open type it need to manually swing the machine body into the horizontal level.






OPEN TYPE

This model allow the operator disassemble and assemble the bowl conveniently, basis on the special Front door design. safety Inspection of the locker and connect with door, ensure the sealing of machine.

C		•	
-S10	ecit	เตลเ	tions
y P	CUII	. Cu .	

Model	Bowl Speed (rpm)	G-Force	Bowl Capacity Nominal (L)	Running Load (kw)	Dimensions (mm)		
					Width	Front-to-Back	Height
105	16300	15700	1200	2.2	600	900	1600
115	16300	16900	1300	3.0	600	900	1600



MIXTURE IN

(4A).

(3) and (4) which exit bowl in separate trays at (3A) and

Typical Application

- Harvesting Biomass •
- Clarify Process Liquid •
- Separating liquid Products Animal
 - Blood Processing
 - Bio remediation •
- Precessing of Granular, Crystalline & Fibrous Suspensions
 - Separation of filterable &
 Non-filterable sludge



Features & Benefits

- Rapid separations
- Consistent Performance
- Range of models facilitate scale-up
- Simple to set up, use and clean for minimal down-time
- No replacement membranes Small
- footprint save space Enclosed
- model available for ambient gas
- control
- For mixtures with very high solid concentration, drum centrifuge can be supplied
- Powerful 3-phase electric motor or customize by various power
- Permanently-lubricated Oversized head bearings maintain long service life

Cooling or Heating Jacket

Tubular separator capable to assemble the cooling or heating jacket, used for some biological material need temperature control.





2017 CATALOGUE

PHILIPPINES

Suite 310 Intramuros Corporate Plaza, Recoletos St., Intramuros, Manila Philippines

Tel: +63 2 2549729

CHINA

No.13093-06-02 Victory Community, Baita District Liaoyang, China

Tel: +86 419 2883215

KOREA

2FL, Hanil bldg, #527-1 Hagi-Dong, Yuseong-Gu, Daejeon, 305-358, Korea

Tel: +82-42-825-4260

AMERICA

6900 46th Street Kenosha, WI 53144, America

Tel:1-262-656-7680